# Extensive chromatin fragmentation improves enrichment of protein binding sites in chromatin immunoprecipitation experiments

Xiaochun Fan, Nathan Lamarre-Vincent, Qian Wang and Kevin Struhl\*

Department of Biological Chemistry and Molecular Pharmacology, Harvard Medical School, Boston, MA 02115, USA

Received June 13, 2008; Revised July 31, 2008; Accepted August 4, 2008

### ABSTRACT

Extensive sonication of formaldehyde-crosslinked chromatin can generate DNA fragments averaging 200 bp in length (range 75-300 bp). Fragmentation is largely random with respect to genomic region and nucleosome position. ChIP experiments employing such extensively fragmented samples show 2- to 4-fold increased enrichment of protein binding sites over control genomic regions, when compared to samples sonicated to a more conventional size range (300-500 bp). The basis of improved fold enrichments is that immunoprecipitation of protein-bound regions is unaffected by fragment size, whereas immunoprecipitation of control genomic regions decreases progressively along with reduced fragment size due to fewer nonspecific binding sites. The use of extensively sonicated samples improves mapping of protein binding sites, and it extends the dynamic range for quantitative measurements of histone density. We show that many yeast promoter regions are virtually devoid of histones.

### INTRODUCTION

Chromatin immunoprecipitation (ChIP) is the standard method to determine the association of proteins with their DNA targets *in vivo* (1). Typically, living cells are treated with formaldehyde to crosslink proteins to DNA, and the crosslinked chromatin is fragmented by sonication. The resulting 'input' sample is immunoprecipitated with an antibody that specifically recognizes the protein of interest, thereby enriching DNA sequences that associate with the protein *in vivo*. The basic measurement of a ChIP experiment is 'IP efficiency' (IE) for a given genomic region, which is defined as the amount of DNA in the IP sample divided by the amount of DNA in the input sample. IP efficiencies for protein-bound regions are higher than for control genomic regions, with the relative level of protein binding being defined as 'fold enrichment' over control regions. For standard ChIP analyses of individual genomic regions, amounts of DNA in the input and IP samples are determined by quantitative PCR, often in real time using an appropriate machine. In large-scale ChIP experiments, the samples are amplified by PCR, and the resulting material analyzed on microarrays (2–4) or by high-throughput DNA sequencing (5–8).

The sensitivity of a ChIP experiment ultimately depends on the ability to separate protein-bound DNA fragments from the background of unbound fragments. In this regard, the quality of the antibody and the IP procedure are critical parameters. However, a different type of background arises from nonspecific association of the protein with genomic DNA that is captured by crosslinking, and such background is unaffected by changes in the IP procedure (9). In principle, the level of nonspecific binding should decrease in accordance with decreasing DNA fragment length, because shorter DNA fragments have fewer nonspecific sites. In contrast, binding to specific target sites should depend primarily on the molarity of the DNA fragments and not the length. Here, we show that extensive fragmentation of chromatin significantly increases fold enrichments of protein binding sites in ChIP experiments.

### MATERIALS AND METHODS

#### **Chromatin preparation**

The yeast experiments utilized a derivative of *Saccharomyces cerevisiae* strain XF236F8 expressing a Pdr1 derivative tagged at the carboxyl-terminus with three copies of the HA epitope. This strain was generated by gene replacement by integrating an appropriate PCR fragment with the *URA3* marker followed by looping out this marker through homologous recombination (10). The procedure for micrococcal nuclease (MNase) digestion of noncrosslinked chromatin was described previously (11).

\*To whom correspondence should be addressed. Tel: +1 617 432 2104; Fax: +1 617 432 2529; Email: kevin@hms.harvard.edu

© 2008 The Author(s)

This is an Open Access article distributed under the terms of the Creative Commons Attribution Non-Commercial License (http://creativecommons.org/licenses/ by-nc/2.0/uk/) which permits unrestricted non-commercial use, distribution, and reproduction in any medium, provided the original work is properly cited.

For ChIP experiments in yeast, cells (40 ml) were grown at 30°C in YPD (1% yeast extract; 2% peptone and 2% dextrose) medium to an optical density of 0.6-0.8 at 600 nm, treated with formaldehyde (1% final concentration) for 20 min at room temperature, and then quenched with 10 ml of 2 M glycine for 5 min at room temperature. The fixed cells were disrupted in a Mini Bead Beater (BioSpec Products, Inc. Bartlesville, OK, USA) with six 3-min cycles at maximum speed, and crosslinked chromatin prepared as described previously (1). The resulting samples were transferred to 2 ml microcentrifuge tubes in 0.5 ml aliquots and sonicated by two methods. First, cells were diluted 2-fold in lysis buffer to a volume of 1 ml and sonicated three times (30 s continuous pulses at about 20% power) in a Branson Sonifier 450 with a microtip probe at  $4^{\circ}$ C, with 2 min cooling on ice between pulses. Second, 0.5 ml cells were sonicated in Misonix Sonicator 3000 (Cole-Parmer Instrument Company, Vernon Hills, IL, USA) with a rotating device and a horn for various times (output level of 6; 10s ON and 10s OFF). Ice and ice-cold water were filled in the chamber to keep samples cool. For both methods, sonicated samples were then centrifuged for 20 min at 16000g to remove insoluble debris. The supernatant were transferred to clean tubes, flash frozen in liquid nitrogen and stored at  $-80^{\circ}$ C.

HeLa cells were grown in MEM supplemented with 10% fetal bovine serum, 100 U/ml penicillin and  $100 \,\mu\text{g/ml}$  streptomycin. Then  $10^7$  cells were trypsinized, fixed in 10 ml MEM supplemented with 1% formaldehyde for 10 min at room temperature, and then quenched by the addition of glycine to a final concentration of 125 mM. After centrifugation, the cells were resuspended in cell lysis buffer (25 mM HEPES pH 7.9, 1.5 mM MgCl<sub>2</sub>, 10 mM KCl, 1 mM DTT, 0.1% NP-40) supplemented with EDTA-free protease inhibitor cocktail (Roche, Indianapolis, IN, USA) and 0.5 mM PMSF and incubated on ice for 10 min. After centrifugation, the crude nuclear pellet was resuspended in nuclear lysis buffer (50 mM HEPES pH 7.9, 140 mM NaCl, 1 mM EDTA, 1% Triton X-100, 0.1% sodium deoxycholate, 0.1% SDS) supplemented with EDTA-free protease inhibitor cocktail and 0.5 mM PMSF. Chromatin was sheared at 4°C, 10 times at 30 second intervals on a Branson Microtip Sonifier 450 (BLD Inc. Garner, NC, USA) set at constant duty and an output level of 4. After centrifugation for 10 min at 16000g, chromatin was sonicated for an additional 0-5 min at 10 s intervals in 0.5 ml aliquots using a cup horn on a Misonix Sonicator 3000 set at level 6.

#### Chromatin immunoprecipitation

Sonicated chromatin preparations were thawed on ice and a portion of the samples saved as input materials. Input materials were decrosslinked at 65°C overnight in the presence of 50  $\mu$ l ChIP elution buffer (50 mM Tris–HCl, pH 7.5; 10 mM EDTA; 1% SDS) and purified through Qiagen columns. DNA was eluted with 30  $\mu$ l TE (10 mM Tris, pH 8.0; 1 mM EDTA), and 10  $\mu$ l of eluate was diluted to 1 ml as input and 20  $\mu$ l of eluate was treated with RNase A/T1 and resolved in 1.5% agarose by electrophoresis. Chromatin from about 0.5 to  $1 \times 10^8$  cells (0.1–0.4 ml chromatin aliquots mentioned above) was incubated with 5µl 8WG16 [for RNA polymerase II (Pol II)] or 10µl F-7 (for HA-tagged protein) for 90 min on a rotating wheel at room temperature, after which 10µl of protein A sepharose was added and incubated for another 90 min. Protein complexes were washed twice with FA lysis buffer (150 mM NaCl), once with ChIP wash buffer (10 mM Tris–HCl, pH 8.0; 250 mM LiCl; 1 mM EDTA; 0.5% Nonidet P-40 and 0.5% sodium deoxycholate) and once with TE. Immunoprecipitated complexes were eluted by incubation with 200µl ChIP elution buffer at 65°C for 1 hour, and samples were decrosslinked at 65°C overnight and purified through QIAgen columns.

#### Quantitation by real-time PCR

For ChIP experiments, input and immunoprecipitated samples were assayed by quantitative PCR to assess the extent of protein occupancy at different genomic regions. PCR reactions contained 2 µl DNA template, 3 µl of  $3.3\,\mu\text{M}$  primer pairs and  $5\,\mu\text{l}$  of  $2\times$  EvaGreen reaction mix (FluoProbes Interchim Montlucon cedex, France). Quantitative PCR was performed on an Applied Biosystem 7500 Fast unit using a 10 min soak at 95°C, followed by 35 cycles of 5s at 95°C, 5s at 53°C and 20 s at 72°C. Threshold cycle (Ct) values were determined at threshold of 0.01. For each amplification product, the IE was determined using the formula  $IE = a * 1.9^{Ct(IP)-Ct(Input)}$ , where a was the constant associated with the ratio and dilution of IP and input materials. Fold occupancy of any given region over an open reading frame (ORF)-free region on chromosome V (chrV) control was determined. IP efficiencies and occupancy values presented here represent an average of three independent experiments and have a standard deviation (SD) of less than  $\pm 25\%$ . Nucleosome scanning analysis was performed as described previously (11). Sequence information of primer pairs used in this study are available upon request.

### RESULTS

# Extensive sonication can reduce chromatin fragments to an average size of 200 bp

Essentially all proteins that associate with specific DNA sequences also associate nonspecifically at a low level with virtually any region of the genome. In the context of ChIP experiments, such nonspecific binding can be seen by the slightly higher IP efficiencies at 'control' regions in IPs with the antibody of interest as compared to control IPs lacking an appropriate antibody (9). As mentioned in the Introduction section, such nonspecific binding should increase with fragment size, whereas specific binding should be largely independent of fragment size. Thus, we reasoned that extensive fragmentation of chromatin should preferentially reduce the nonspecific binding, and hence might result in improved fold enrichments in ChIP experiments.

Our heretofore standard procedure for yeast cells utilized a Branson Sonifier 450 with a microtip probe for 90 s,



**Figure 1.** Comparison of yeast chromatin fragment size with different sonication time and methods. (A) Chromatin samples prepared by sonication in Misonix Sonicator 3000 with a horn or in Branson Sonifier 450 with a microtip probe were removed at various time, decrosslinked, purified, treated with RNase and resolved in 2% agarose gel. Lanes with molecular weight markers are labeled M1 (PCR marker) and M2 (100 bp DNA, both from New England Biolabs Ipswich, MA, USA) and lanes with sonicated chromatin are labeled with length of time (in minutes) or 90 s with a microtip probe (P). Sizes of molecular weight markers are indicated. (B) DNA prepared from crosslinked chromatin sonicated for 16 min with a horn (Chr) or from mono-nucleosomes prepared by MNase treatment of noncrosslinked chromatin (MN) were analyzed with the indicated primers pairs spanning the *PHO5* gene. Values are normalized to those obtained with genomic DNA and presented with respect to the average value (defined as 1) in the same sample.

resulting in DNA fragments in the range of 300–500 bp, which is typical for ChIP experiments described by other laboratories. We were able to reduce DNA fragment size further using a Misonix Sonicator 3000 with a horn (Figure 1A). Sonication for 4 min yielded a size range of DNA fragments similar to that achieved with the Branson Sonifier 450, but sonication for 16 min yielded DNA fragments ranging from 75 to 300 bp in length. In addition, the yield of solubilized chromatin increased about 2- to 3-fold in samples sonicated for longer time (data not shown). Sonication times as long as 60 min did not further reduce DNA fragment size significantly.

# Fragmentation is largely random with respect to genomic region and nucleosome position

The average size of the extensively sonicated chromatin is roughly similar to that of mono-nucleosomes, and we considered the possibility that fragmentation would preferentially occur within the relatively flexible linker regions between nucleosomes. However, nucleosome-scanning analysis (11) across the *PHO5* promoter and coding region showed that the extensively sonicated sample has constant level of DNA at all positions tested. In contrast, parallel analysis of mono-nucleosomal DNA generated by MNase treatment of noncrosslinked chromatin from yeast cells showed the expected pattern of positioned nucleosomes (Figure 1B). Thus, extensive sonication of crosslinked chromatin generates fragmentation patterns that are largely random, although subtle, locus-specific preferences cannot be excluded.

## Increased fold enrichments upon increased fragmentation of chromatin

We performed ChIP analysis for RNA Pol II on samples generated at various times of sonication. As shown in Figure 2A, for genomic regions tested that show high Pol II occupancy, the various samples showed similar IP efficiencies. In contrast, at control genomic regions, the IP efficiency decreased significantly in samples with longer sonication time. When IP efficiencies were normalized to a negative control, an ORF-free region on chrV, fold enrichments at all Pol II-associated regions increased progressively in accordance with the time and extent of sonication, whereas fold enrichment at the telomeric region (TEL) was unaffected (Figure 2B). As expected from the similar distributions of DNA fragment sizes, the ChIP results obtained from samples prepared with the Branson Sonifier 450 were comparable to those of samples prepared sonication for 4 min with the Misonix Sonicator 3000. Thus, by reducing DNA fragment sites to 75–300 bp (16 min sonication with the horn) from the more conventional 300–500 bp (Branson Sonifier 450), we were able to obtain fold enrichments that are about 2-fold higher than that obtained previously.

# PCR product length does not significantly affect fold enrichments

One issue with samples containing reduced size of chromatin fragments is that quantitative analyses involving primer pairs producing long PCR products might be adversely affected. In principle, the input and immunoprecipated samples should be equally affected by the length of the PCR products used in the analysis, thereby having no effect on fold enrichments. To test this, we compared the fold enrichments at *POL1* and *BTN2* using primer pairs that generated small (74 and 64 bp) or large (219 and 349 bp) PCR products (Table 1). In terms of fold enrichment, there was no significant difference between each pair of primer set in chromatin samples prepared with various sonication time and different methods (Figure 2C).



Figure 2. RNA Pol II occupancy at some genomic regions. Pol II occupancy at some genomic regions for chromatin preparations with different sonication time and methods was presented as IP efficiency (A) or fold over ORF-free region on chrV (B). (C) Pol II occupancy determined for genes with primer pairs giving different sizes of PCR products. The size of PCR products can be found in Table 1. (D) Differences of Ct values between primer pairs giving different sizes of PCR products. Error bars show SD from three independent experiments.

**Table 1.** Oligos used to compare the size effect of PCR products on fold enrichments in Figure 2D and to map Pdr1 binding sites at *VHR1* promoter in Figure 3C, with the coordinates relative to ATG

Primer	5′	3′	Size
POL1(S)	+713	+786	74
POL1(I)	+2499	+2717	219
BTN2 (S)	+657	+720	64
BTN2 (I)	+559	+907	349
Α	-1010	-943	68
В	-962	-891	72
С	-881	-812	70
D	-752	-690	63
Е	-657	-597	61
F	-558	-499	60
G	-484	-425	60
Н	-345	-286	60

As expected, however, it took more PCR cycles (higher Ct values) to reach the measurement threshold in samples that were more extensively sonicated (Figure 2D), reflecting the decreasing population of DNA molecules that are long enough to include sequences corresponding to both primers.

# Improved mapping of protein binding sites upon extensive fragmentation of chromatin

A protein bound to a specific DNA site *in vivo* generates a predicted ChIP profile that depends on the size of the fragmented DNA in the input chromatin sample and the length of the PCR products used for the analysis (12). Peak ChIP signals are centered at binding sites and extend the length of the PCR product, whereupon gradually and symmetrically decrease at regions flanking the peak in a manner that depends on DNA fragment length. Thus, chromatin samples with shorter DNA

fragments should yield ChIP profiles with sharper peaks, and hence improve the resolution of the binding site.

To examine the effect of fragment size on mapping protein binding sites, we examined the association of (HA)<sub>3</sub>tagged Pdr1. In accordance with the Pol II results above, IP efficiencies of regions containing Pdr1 binding sites (PDR5, SNQ1 and PDR16 promoters) were unaffected by sonication times, whereas they significantly decreased at control genomic regions (Figure 3A). As a consequence, fold enrichments increased as a function of sonication time, with the most extensively fragmented samples showing 4-fold higher ChIP signals than that obtained with the conventional method using a sonication microtip probe (Figure 3B). We mapped Pdr1 occupancy at the VHR1 promoter region using eight PCR primer pairs spanning a 700 bp region encompassing the VHR1 promoter (Table 1). For chromatin samples prepared by sonication for 16 min with a horn, Pdr1 association appeared as a sharp peak at -752 to -690; while Pdr1 profile was significantly flatter for chromatin samples prepared by sonication for 90s with a microtip probe (Figure 3C). These results are consistent with the putative Pdr1 binding sites near -736 and -728 at the VHR1 promoter (13). Similar results were observed at the PDR5, PDR16 and TPO1 promoters (data not shown).

#### Improved measurements of histone density upon extensive fragmentation of chromatin reveals that many yeast promoter regions are virtually devoid of histones

Although most DNA in eukaryotic cells exists in the form of nucleosomes, many promoter (and perhaps other) regions appear to be depleted of histones (11,14–17). However, it has been difficult to quantitate the extent of nucleosome depletion. First, analyses involving mononucleosomal DNA generated by MNase treatment are complicated by the strong DNA sequence preferences of



Figure 3. Pdr1 occupancy at some genomic regions. Pdr1 occupancy at some genomic regions for chromatin preparations with various sonication time and different methods were presented as IP efficiency (A) or fold over chrV (B). (C) Pdr1 occupancy across *VHR1* promoter, presented as fold over chrV, with values normalized to the maximum value (defined as 100) in the same sample. Error bars show SD from three independent experiments.

MNase (18) that preferentially cleave AT-rich promoter regions. Second, measurements of histone density by ChIP using antibodies against histones (or epitope-tagged histones) are problematic because the 300–500 bp DNA fragments from standard ChIP samples will have 2–4 nucleosomes and histone crosslinking to DNA is extremely efficient. In this regard, the 2-fold reduction in histone density at many yeast promoter regions (11) may be a significant underestimate of nucleosome depletion, especially when nucleosome-free regions are shorter than the size of the chromatin fragments.

To address this issue, we performed a ChIP experiment on the extensively sonicated chromatin sample using an antibody against histone H3. As shown in Figure 4, many yeast promoter regions are virtually devoid of histones, with histone densities <10% of those observed at coding regions. As expected, the *PHO5* (19) and *GAL1,10* (20) promoter regions have much higher histone densities, although these are typically 2-fold below that of coding regions. Also, as expected (16,21), coding regions of highly active genes (*RPL3* and *ENO2*) also show approximately 2-fold lower levels of histone density in comparison to those observed at a poorly transcribed gene (*POL1*), ORF-free region and TEL. Thus, the use of extensively sonicated samples greatly increases the



Figure 4. Histone density at yeast promoter and coding regions. Histone H3 occupancy at the indicated promoter and coding regions, presented as IP efficiency. Error bars show the SD from three independent experiments.

quantitative accuracy and dynamic range of histone density experiments.

# Improved fold enrichment upon extensive fragmentation of crosslinked chromatin from human cells

To extend these results to ChIP experiments in mammalian cells, crosslinked chromatin from HeLa cells was prepared by the standard procedure using the Branson Sonifier 450. This material was then sonicated for an additional 2 or 5 min with the Misonix Sonicator 3000, resulting in a progressive decrease in the size of DNA (Figure 5A). After 5 min, the bulk of the DNA was between 150 and 250 bp, which is slightly larger than the size range obtained for yeast cells. For the experiment shown, generating DNA in the 150–250 bp range required this two-step sonication procedure; treatment with the Misonix Sonicator 3000 alone was insufficient. However, the precise conditions will likely vary among cell types.

We performed ChIP analysis for RNA Pol II on the above samples with the same antibody used in the yeast experiments. In accordance with the results in yeast, longer sonication times resulted in more dramatic reduction in IP efficiencies at control genomic than at the *ING1* promoter and *EEF1A1* coding regions (Figure 5B). As a consequence, enrichment at the *ING1* promoter and *EEF1A1* coding region was 3- to 4-fold higher in the most extensively sonicated sample as compared to the sample generated with the Branson Sonifer 450 alone (Figure 5C). Attempts to improve the fold enrichment by more extensive sonication (10 min with the Misonix Sonicator 3000) resulted in reduced IP efficiency and fold enrichment at the *ING1* promoter, and we suspect that this is due to damage to Pol II.

#### DISCUSSION

Extensive sonication of formaldehyde-crosslinked chromatin samples from both yeast and human cells can reduce average fragment size to  $\sim 200$  bp (range 75–300 bp). Such extensively fragmented chromatin samples yield a 2- to 4-fold increase in fold enrichments in ChIP experiments when compared to samples sonicated to a more conventional size range (300–500 bp). The basis of improved fold enrichments is that IP efficiencies at protein binding sites are unaffected by fragment size,



Figure 5. RNA Pol II occupancy in human cells. (A) Chromatin samples prepared by sonication in Branson Sonifier 450 with a microtip probe and subsequent sonication with the Misonix Sonicator 3000 for the indicated time. After decrosslinking, the resulting DNAs were resolved in a 2% agarose gel along with molecular weight markers. (B) IP efficiencies at the *ING1* promoter, *EEF1A1* coding and control (Ctr) genomic regions for the indicated samples. (C) Fold enrichment of Pol II occupancy at the *ING1* promoter and *EEF1A1* coding regions relative to the average of the control genomic regions.

whereas IP efficiencies at control genomic regions decrease progressively along with reduced fragment size. These experimental observations are in accord with the prediction that the level of nonspecific binding should decrease in accord with decreasing DNA fragment length, because shorter DNA fragments have fewer nonspecific sites, whereas binding to specific target sites should depend primarily on the molarity of the DNA fragments. We do not know whether fragment size affects the background in ChIP experiments due to the population of molecules that are unbound by the protein of interest.

The use of extensively fragmented chromatin samples should increase the sensitivity of all ChIP applications, because sensitivity is essentially defined by fold enrichment. First, extensively fragmented samples will be particularly beneficial for ChIP experiments involving proteins that inefficiently crosslink to chromatin. For example, ChIP experiments that give rise to 3- to 5-fold enrichments by conventional procedures would be expected to yield 6to 20-fold enrichments, a very significant difference. Second, increased sensitivity of the ChIP assay will also be beneficial when proteins associate weakly with genomic regions, either due to an intrinsically weak binding site or as a consequence of environmental or genetic conditions. Third, highly fragmented chromatin samples improve resolution of where the protein associates with DNA. Fourth, extensive sonication greatly improves quantitative measurements of histone density, and it reveals that many yeast promoter regions are virtually devoid of histones. Lastly, extensively fragmented chromatin samples should significantly improve the ability to distinguish target versus nontarget sites in ChIP-chip and ChIPsequencing experiments. In both of these large-scale applications, variability due to the amplification procedure, cross-hybridization and counting error due to the limited amount of DNA sequencing decreases the sensitivity of the assay. The higher chromatin yields arising from extensive sonication will help reduce variability due to amplification and counting error. In addition, extensive sonication produces DNA fragments that are of ideal size for making libraries for ChIP sequencing.

Unlike chemical or enzymatic treatments, sonication is an idiosyncratic procedure that depends on many variables and needs to be calibrated empirically. In our laboratory, we were able to fragment yeast chromatin to a size range of 75-300 bp using either the Misonix Sonicator 3000 with a horn for 16 min or the Branson Sonifier 450 with a microtip probe for 8 min. One advantage of the Misonix Sonicator is the ability to process up to 20 samples simultaneously, and the use of a rotating device minimizes variation in fragment size among samples. The degree of sonication might also be affected by crosslinking time, formaldehyde concentration and the concentration of detergent in the broken cells, although we have not investigated these parameters systematically. In our experiments with the HeLa cell derivative, fragmentation to the 150-250 bp range required a two-step sonication procedure involving both the Branson Sonifier 450 and the Misonix Sonicator 3000. However, as chromatin fragmentation by sonication can differ among mammalian cell types, the procedure needs to be optimized for each cell type. In addition, we found that excessive sonication reduced IP efficiency and fold enrichments, probably due to damage to the protein of interest.

In summary, our results indicate that chromatin can be sonicated to a size range of 75–300 bp, and that such extensive sonication is highly desirable for ChIP experiments. Thus, we believe that it is well worth the effort to optimize sonication conditions that permit such extensive fragmentation of chromatin.

### ACKNOWLEDGEMENTS

We thank Joseph Fleming and Marianne Lindahl Allen for communicating unpublished results, Heather Hirsch and Zarmik Moqtaderi for technical advice and Joseph Geisberg for useful discussions.

### FUNDING

Damon Runyon-Walter Winchell Cancer Research Foundation postdoctoral fellowship (Q.W.); National Institutes of Health (GM30186 to K.S.). Funding for open access charge: NIH GM 30186, GM 53720, and CA 107486.

Conflict of interest statement. None declared.

### REFERENCES

- Aparicio,O.M., Geisberg,J.V., Sekinger,E.A., Yang,A., Moqtaderi,Z. and Struhl,K. (2005) Chromatin immunoprecipitation for determining the association of proteins with specific genomic sequences in vivo. *Curr. Protoc. Mol. Biol.*, 69, 21.23.21–21.23.17.
- Bernstein, B.E., Humphrey, E.L., Liu, C.L. and Schreiber, S.L. (2004) The use of chromatin immunoprecipitation assays in genomewide analyses of histone modifications. *Methods Enzymol.*, 376, 349–360.
- Buck, M.J. and Lieb, J.D. (2004) ChIP-chip: considerations for the design, analysis, and application of genome-wide chromatin immunoprecipitation experiments. *Genomics*, 83, 349–360.
- Ren,B. and Dynlacht,B.D. (2004) Use of chromatin immunoprecipitation assays in genome-wide location analysis of mammalian transcription factors. *Methods Enzymol.*, 376, 304–315.
- Impey,S., McCorkle,S.R., Cha-Molstad,H., Dwyer,J.M., Yochum,G.S., Boss,J.M., McWeeney,S., Dunn,J.J., Mandel,G. and Goodman,R.H. (2004) Defining the CREB regulon: a genome-wide analysis of transcription factor regulatory regions. *Cell*, **119**, 1041–1054.
- 6. Johnson, D.S., Mortazavi, A., Myers, R.M. and Wold, B. (2007) Genome-wide mapping of in vivo protein-DNA interactions. *Science*, **316**, 1497–1502.

- 7. Wei,C.L., Wu,Q., Vega,V., Chiu,K.P., Ng,P., Zhang,T., Shahab,A., Yong,H.C., Fu,Y.T., Weng,Z. *et al.* (2006) A global map of p53 transcription factor binding sites in the human genome. *Cell*, **124**, 207–219.
- Robertson,G., Hirst,M., Bainbridge,M., Bilenky,M., Zhao,Y., Zeng,T., Euskirchen,G., Bernier,B., Varhol,R., Delaney,A. *et al.* (2007) Genome-wide profiles of STAT1 DNA association using chromatin immunoprecipitation and massively parallel sequencing. *Nat. Methods*, 4, 651–657.
- Struhl,K. (2007) In Zuk,D. (ed.), Evaluating Techniques in Biochemical Research. Cell Press, Cambridge, MA, pp. 29–33.
- 10. Schneider, B.L., Seufert, W., Steiner, B., Yang, Q.H. and Futcher, A.B. (1995) Use of polymerase chain reaction epitope tagging for protein tagging in *Saccharomyces cerevisiae*. *Yeast*, **11**, 1265–1274.
- Sekinger, E.A., Moqtaderi, Z. and Struhl, K. (2005) Intrinsic histone-DNA interactions and low nucleosome density are important for preferential accessibility of promoter regions in yeast. *Mol. Cell*, 18, 735–748.
- Kadosh,D. and Struhl,K. (1998) Targeted recruitment of the Sin3-Rpd3 histone deacetylase complex generates a highly localized domain of repressed chromatin *in vivo*. *Mol. Cell. Biol.*, 18, 5121–5127.
- Teixeira, M.C., Monteiro, P., Jain, P., Tenreiro, S., Fernandes, A.R., Mira, N.P., Alenquer, M., Freitas, A.T., Oliveira, A.L. and Sa-Correia, I. (2006) The YEASTRACT database: a tool for the analysis of transcription regulatory associations in *Saccharomyces cerevisiae*. *Nucleic Acids Res.*, 34, D446–D451.
- 14. Ng,H.H., Robert,F., Young,R.A. and Struhl,K. (2003) Targeted recruitment of Set1 histone methylase by elongating Pol II provides a localized mark and memory of recent transcriptional activity. *Mol. Cell*, **11**, 709–719.
- Bernstein, B.E., Liu, C.L., Humphrey, E.L., Perlstein, E.O. and Schreiber, S.L. (2004) Global nucleosome occupancy in yeast. *Genome Biol.*, 5, R62.
- Lee,C.K., Shibata,Y., Rao,B., Strahl,B.D. and Lieb,J.D. (2004) Evidence for nucleosome depletion at active regulatory regions genome-wide. *Nat. Genet.*, 36, 900–905.
- Yuan,G.-C., Liu,Y.-J., Dion,M.F., Slack,M.D., Wu,L.F., Altschuler,S.J. and Rando,O.J. (2005) Genome-scale identification of nucleosome positions in *S. cerevisiae. Science*, **309**, 626–630.
- 18. Horz, W. and Altenburger, W. (1981) Sequence specific cleavage of by micrococcal nuclease. *Nucleic Acids Res.*, 9, 2643–2658.
- Almer, A., Rudolph, H., Hinnen, A. and Horz, W. (1986) Removal of positioned nucleosomes from the yeast *PHO5* promoter upon *PHO5* induction releases additional upstream activating DNA elements. *EMBO J.*, 5, 2689–2696.
- Cavalli,G. and Thoma,F. (1993) Chromatin transitions during activation and repression of galatose-regulated genes in yeast. *EMBO J.*, **12**, 4603–4613.
- Schwabish, M.A. and Struhl, K. (2004) Evidence for eviction and rapid deposition of histones upon transcriptional elongation by RNA polymerase II. *Mol. Cell. Biol.*, 24, 10111–10117.